

DW-GI Series **Vertical Autoclave**

Advanced

- **Colorful LCD screen**

LCD screen with English and Chinese choice, using text, data, and graphical methods to dynamically display various statuses, the sterilization process is clear at a glance.

- **6 level of exhaust speed control (TW, TR only)**

After sterilization, you can choose different exhaust speed, and you can adjust it freely at any time of exhausting.

- **Condenser (Option)**

Equipped with a professional condenser, the steam generated during the sterilization process is converted into water through the condenser, achieving zero discharge of water vapor/mist.

- **F0 value calculation**

F0 value is automatically calculated, can print out by micro-printer.

- **Drying (TR only)**

Fully automatic drying: Equipped with an air pump and a 0.2μm air filter, it automatically starts drying after sterilization.

After sterilization, clean air will be put in to replace the wet air inside, so as to improve drying efficiency.

(Fully automatic drying is an option for TR only)

- **Biosafety function (TWS, TRS only)**

Use a 0.2 μm PTFE filter packed in a housing to ensure that no biohazard aerosols, microorganisms and viruses are passing out through the filter during exhausting, so as to make sure the biology safety of the lab.

- **5 level access permission**

The operation permissions for instruments are divided into five levels, with each level of user having an independent password and multiple IDs that can be assigned. Each level of user has limited access to certain usage permissions.

- **Fast cooling fan**

To shorten cooling process after sterilization.

- **Auto water fill (Option)**

Offer pipe and pump to connect to the direct water supply to achieve automatic water feeding, must be added while production.

- **Air inlet filter (TR as option)**

Filtration is higher than 99.99% in particle with diameter of 0.2 μm.

- **Monitor of exhaust filter (TWS, TRS only)**

Exhaust filter is put inside a pressure resistant housing (stainless steel material) with an extra PT100 temperature sensor to monitor the temperature, so as to make sure the filter and housing are also being sterilized thoroughly, and will not lead to secondary biohazard pollution.

- **Fully automatic drying (TR only)**

Equipped with an air pump and a 0.2μm air filter, it automatically starts drying after sterilization. After sterilization, clean air will be put in to replace the wet air inside,

Intelligent

1. Calibration port

Offer temp.and pressure interface, together with the 3Q adapter joints, can insert max 15 sensors one time to do the calibration work.

2. Program memory storage system

User can set their own program, and saved in

3. Saturated steam monitoring

The system automatically monitors the purge out of cold air to ensure the pure steam environment.

4. Special software for drug packing materials(option)

The product meets the 2020 Chinese Pharmacopoeia or EP10.2 standards.

Safety Features

- **Nice look:**The lid and bench corners are covered with insulation to avoid getting scald

Inter lock:Only when the temperature and pressure are safe inside, can the door be opened. If the door is not closed in place, it cannot start working.

- **Double protection for over pressure:**Mechanical pressure gauge and electric sensing system, once pressure is abnormal, system will cut power, release the pressure and send out alarm.
- **Safety test function:**Regular testing of system security and safety protection measures is possible.
- **Electrical production:**Prevention for over current,short ciucuit and electricity leakage.
- **Lid closing check:**Machine can not start if lid is not locking well by checking mechanically and electrically.
- **Triple protection for water shortage:**The bottom of the sterilization chamber is equipped with three different dry burning protection devices: liquid expansion type, copper temperature sensing type, and electric ion type(water level sensor), to avoid misjud-gment caused by a single method.
- **Auto trouble shooting:**System will monitor the working status, and will cut the power, send out alarm, and report the error information once any abnormal situation detected.
- **Temperature monitor:**Power cut and alaraming sent, if there is temperature is too high or temperature lising is abnormal.
- **Saturated steam monitoring:**The system automatically monitors the purge out of cold air to ensure the pure steam environment.

DW-GITW models
DW-GI29TW/
DW-GI54TW/
DW-GI80TW/
DW-GI100TW



DW-GITR models
DW-GI54TR/
DW-GI80TR/
DW-GI100TR

Hommization

1. Just right depth

Conforming to ergonomics, the chamber depths of 29/54/80/100Lchamber are 378/725/670/780 (mm), making it easy to take and place samples and maintain the instrument

2. Front water tank

Front water tank, easy to plug in and out. (TW only)

Specifications

Model \ Function	DW-GITW	DW-GITR
6 level exhaust speed control	Standard	Standard
Auto water fill	Option	Option
Condensor	Option	Option
Dedorization		Standard
Drying	——	Standard(Semi-automatic drying); Option(fully automatic drying)
Cooling fan	Standard	Standard
Access permission	5 level	5 level
USB port	Standard	Standard
Sterilization record	Standard	Standard
Printing set	Option	Option
Micro-printer	Option	Option
Display system	LCD	LCD
Load thermomete(floating PT100 sensor)	Option	Option
Auto drainage	Option	Option
Steam collecting way	Internal	Internal
Prevacuum		——
Auto cleaning	Option	Option
Digital pressure sensor	Option	Option
Inlet air filter	——	Option
Air pump	——	Option
S Biosafety function	DW-GI-TWS only	DW-GI-TRS only

Note: Requirements for sterilization water: distilled water (conductivity not less than 3μS/cm)
When configured as an S biosafety type, a condenser is standard

Parameter

Model		Standard	Standard	DW-GI54TW	DW-GI80TW	DW-GI100TW
		Drying	DW-GI29TW	DW-GI54TR	DW-GI80TR	DW-GI100TR
Difference			DW-GI-TW:Fully automatic with cooling fan DW-GI-TR:DW-GI-TW+semi automatic drying, if need fully automatic drying need to add 530USD			
Capacity(L)			29	54	80	100
Dimension (LxWxH)(mm)	TW		460x542x976	460x542x1105	620x745x1090	620x745x1200
	TR			460x542x1105	620x726x1090	620x730x1200
Chamber dimension (Dia*H)(mm)			Φ320x378	Φ320x725	Φ400x700	Φ400x800
Rated power(W)			2300	2600	4600(2900 as option)	4600
Sus304 basket size(pcs)			(Φ298x260)x1	(Φ298x195)x3	(Φ360x275)x2	(Φ360x215)x3
Power requirement			220V±10% 16A 50HZ/60HZ		220V±10% 32A 50HZ/60HZ	
Temperature diTPlay accuracy			0.1℃			
Working environment			5℃~40℃ Relative humidity ≤85%			
Chamber material			SUS304 stainless steel (316L as option)			
Sterilization temp.range			105℃~135℃			
Sterilization time range			1~6000 min			
Melting temp. range			60℃~115℃			
Melting time			1~6000 min			
Warming temp.			45℃~79℃			
Warming time			1~9999 min			
Memory function			60 program saved in system,999 cyeles saved by USB			
Appointment timer range			Auto start up appointment)0～15 days			
Exhaust TPeed control			Fully automatic internal discharging with 6 levels choice			
Sterilization mode			Liquid,liquid with warming,solid,wrapped instrument,fabric,rubber, fast,waste,agar,self-defined Drying(TR only)			
Drying(TR、 TP only)			Drying time:1~300min Drying temp.:80~160℃			
Controller			"SMART III" smart microcomputer controller			
Cooling way			Fast cooling fan			
Pressure gauge range			-0.1~0.5MPa			
Pressure valve set pressure			0.29Mpa			
Vacuum pump			Pulse times: 1~9 times Vacuum degree:≤-80kpa Vacuum leakage rate (per min)<0.13kpa			
Safety device			Lid closing checking,dry scorch protection,,safety valve,water level sensor, over temperature/pressure protection,temperature/pressure monitor, electricity safety(over current/short circuit/leakage),cooling lock,anti-scald lid and bench corner,automatic troubleshooting			

Note :Baskets can also be customized, as well as melting and warming temperature range.

Spare part option




Micro printer
To print temp.,time and F0 value



Digital pressure sensor



Load thermometer
Floating PT100 sensor



Tray



Stainless steel bucket



IQ/OQ/PQ
(3Q)adapter joint



Adjust foot

Printing set:printer+digital
pressure sensor:can print temp. and
pressure in data and flow chart way

IQ/OQ/PQ documents