

Tabletop Low-Speed Large Capacity Centrifuge TD6

TD6 is low speed Large Capacity centrifuge, Which is the idea machine for Biological, hospital clinics laboratories, Epidemic prevention stations and to separate samples of blood, plasma, Urine and so on.

Features:

- ♦ Stainless Steel Chamber, Full metal structure and microprocessor controlled
- Large capacity and many Rotors for choose, to meet the requirement of user
- ◆ Safety Interlock system that prevents opening the centrifuge lid during running.
- With air cooling fans to avoid high temperature with long time working.
- Digital or LCD display for choose, which indicated the Speed, time, RCF etc.
- ◆ 10 programs can be stored, quickly recall the last used program.
- ♦ Imbalance sensor, Over speed protection
- With Emergency lid lock release, in case of power off.













Parameters

Model	TD6		
Max Speed	6000r/min		
Speed Accuracy	± 30 rpm		
Max RCF	5600×g		
Timer Range	0~99min		
Noise	≤65dB(A)		
Power Supply	AC220V & 110V 50Hz 5A		
Dimension	$602\times502\times412\mathrm{mm}(\mathrm{L}\times\mathrm{W}\times\mathrm{H})$		
Weight	52kg		
Rotor No.	Rotor Capacity	Speed (r/min)	RCF (xg)
No.1 Angle Rotor	12×10/15ml	5000	4390
No.2 Swing Rotor	24×5ml	4200	3120
No.3 Swing Rotor	108×5ml	4000	4200
No.4 Swing Rotor	16×10/15ml	5000	4100
No.5 Swing Rotor	8×50ml	4200	2120
	4×50ml		3120
No.6 Swing Rotor	8×100ml	4200	3120
No.7 Swing Rotor	4×500ml	4000	3600
No.8 Well Plates Rotor	2×2×96holes	4000	2300