

Micro Hematocrit Centrifuge TG12-W

Features:

• TG12-W capillary blood centrifuge is mainly used to determine the value of the blood hair cells, as well as the centrifuge of the

micro blood and micro capacity. Full metal structure and bushless motor.





Parameters

Model



Max speed	12000r/min
Max RCF	15800xg
Rotor Capacity	24 Pieces of Capillary
Timer Range	0-99min
Motor	Microprocessor control, DC Brushless motor
Noise	≤55dB(A)
Power Supply	AC220V&110V 50Hz 2A
Dimension	380*280*200mm(L*W*H)
Weight	15.5KG

