## DRAWELL Artist of Science

## 0.01g TD-A LCD Electronic Balance (External Calibration)



High precision loadcell, compact design.

Full plastic housing.

Stainless steel Plstter.

Super bright LCD display with backlight.

Mains adapter supplied as standard.

Height adjustable feet.

Full Capacity substraction.

Simple 5 key panel for easy operation.

Stabilization time 1.5 seconds.

Overload Protection.

Below Balance Weighing Facility.

Selectable measure units, g,oz,ct.



Checking Weighing, Piece counting function

## Specifications

Model	<b>TD1002A</b>	<b>TD2002A</b>	<b>TD3002A</b>	<b>TD5002A</b>		
Capacity	<b>0-100g</b>	<b>0-200g</b>	<b>0-300g</b>	<b>0-500g</b>		
Readability	0.01g					
Repeatability	±0.02g					
Linearity	±0.02g					
Net Weight	1200g					
<b>Power Supply</b>	DC 6V-9V (AAA battery*4pcs)					

Model	<b>TD10002A</b>	<b>TD20002A</b>	<b>TD30002A</b>	<b>TD50002A</b>		
Capacity	<b>0-1000g</b>	<b>0-2000g</b>	<b>0-3000g</b>	<b>0-5000g</b>		
Readability	0.01g					
Repeatability	±0.02g					
Linearity	±0.02g					
Net Weight	1200g					
<b>Power Supply</b>	DC 6V-9V (AAA battery*4pcs)					



www.drawell.com.cn