

0.1g TD-A LCD Electronic Balance (External Calibration)

Features

- High precision loadcell, compact design.
- Full plastic housing.
- Stainless steel Plstter.
- Super bright LCD display with backlight.
- Mains adapter supplied as standard.
- Height adjustable feet.
- Full Capacity substraction.
- Simple 5 key panel for easy operation.
- Stabilization time 1.5 seconds.
- Overload Protection.
- Below Balance Weighing Facility.
- Selectable measure units, g,oz,ct.
 Checking Weighing, Piece counting function.



Specifications

Model	TD2001A	TD6001A	TD-10001A	TD-20001A	TD-30001A	TD50001A
Capacity:	0-200g	0-600g	0-1000g	0-2000g	0-3000g	0-5000g
Readability:	0.1g	0.1g	0.1g	0.1g	0.1g	0.1g
Repeatability:	±0.2g	±0.02g	±0.2g	±0.2g	±0.2g	±0.2g
Linearity:	±0.2g	±0.02g	±0.2g	±0.2g	±0.2g	±0.2g
Platter Size:	155*143mm					
Dimensions:	280mm x 180mm x80mm					
Net Weight:	1200g					
Power Supply:	DC 6V-9V(AAA battery *4pcs)					