

UV analyzer DW-UVA-7 (UV Darkroom)

Description

◆ DW-UVA-7 UV analyzer is used in many fields of nucleic acid examination, including molecular biology, biochemistry, medicine and forensic examination, as well as biological agents, such as observation of fluorescence.

Features

- ♦ Instrument is installed 254 nm and 365 nm uv light source,UV light bulb be arranged closely.
- Fully enclosed with high quality Darkroom, No need to operate in darkroom. Convenient to use, simple to operate.
- Equipment can be installed digital camera to shoot and take photos.
- The built-in reflector enhances UV illumination.





Specification

Transmission wavelength	254nm、365nm
Transmission Filter Size	200mmX50mm
Instrument dimension	280×70×120mm
Instrument weight	3.2kg
Power	DC 220V,50Hz