

## UV Analyzer DW-UVA-5 (UV lamps)

## Description

◆ UVA-5 UV analyzer are used in many fields of nucleic acid examination, including molecular biology, biochemistry, medicine and forensic examination, as well as biological agents, such as observation of fluorescence.

## Features

- ♦ The design is simple and smooth, easy to operate and portable, ergonomic, switch is very convenient.
- ♦ With a single and double wave two series can be selected, 365nm and 254nm or double 302nm combination
- ♦ The built-in reflector enhances UV illumination.



## Specification

Transmission wavelength	254nm、365nm OR double 302nm
Transmission Filter Size	200mmX50mm
Instrument dimension	280×70×120mm
Instrument weight	1.2kg
Power	DC 220V,50Hz