

## DTF-8D Water Bath

## Features

- Stainless steel chamber and top cover.
- ♦ 3 independent water tanks and temperature controls. Can separately set and control temperature.
- Freely set operating temperature, with alarm function.
- LCD display for set temperature and measured temperature.
- Equipped with leakage protection.



## Specifications

Model	Power Supply	Temperature Range (°C)	Power Rating (W)	Temperature Fluctuation (°C)	Chamber Size (W×D×H)cm	Exterior Size (W×D×H)cm	Package Size (W×D×H)cm	Net/Gross Weight (kg)
DTF-8D	220V, 50HZ	RT+5~99	800	±0.5	(14×16.8×11.5)×3	50×30×31	58×35×38	9/14