



DW-T3C Binocular Stereo Microscope

- ◆ This series microscope is a basic two-steps fixed stereo microscope. It is with large clear range, uniform field and easy operation. This microscope is used extensively world agriculture, as well as in electronic component manufacturing, specially suited for LED, PCB inspection, pressing plating and electronic component inspection.

Eyepiece	Wide Field Plan-scope Eyepiece: WF10X/18mm
Objective	Grouped From Any Two Objectives of 1X, 2X, 3X and 4X, Working Distance: 1X-55mm, 2X-81mm, 3X-61mm, 4X-57mm
Head	Articulated Free Binocular Head, <u>Interpupillary</u> Distance: 55-75mm, 45°Inclined
Focusing Adjustment	Rack and Pinion Focusing with Mechanism for Adjusting Focusing Torque
Arm Type	Pole
Illumination	Top and Bottom light, Light Optional: 12V/10W halogen lamp, 1W LED
Optional	1.WF15X, WF20X; 2. Rechargeable batteries
Net Weight	2.4kg

