

LCD Digital 3D Shaker

Model: SK-D3309-Pro

3D shaker is used for uniform shaking mixing in flasks, culture flasks, petri dishes and tubes. ideal for culturing a variety of cell cultures, DNA extractions, distributing of cells and gel staining & destaining of cells.

Features

- Speed range of 10-70 rpm
- Maximum load capacity of 5kg with platform
- Independent LCD display of time and speed
- Incubator compatible
- Various platform accessories
- Three-dimensional rocking motion with an angle of 9°



Specifications

Model#	SK-D3309-Pro
Tilt angle	9°
Max. load capacity (with platform)	5kg
Motor type	DC motor
Motor input power	40W
Motor output power	24W
Speed range	10-70rpm
Speed display	LCD
Timer display	LCD
Time settings range	1min-19h59min
Voltage	100-240V, 50/60Hz
Power	50W
Weight	9kg
Permissible ambient temperature and humidity	5-40°C 80%RH
Dimension[WxDxH]	360×430×106mm
Protection class	IP21