

## Electric Heating Double Distillation

### DZ5C/DZ10C/DZ20C

#### Performance and characteristic:

- All of products are made from stainless steel.
- Have two container, and firm reasonable.
- Condenser is made of seamless stainless-steel pipes and so with very good cooling effects & great capacity of water produce.
- When lack water, will auto stop heating, when the return water reach the normal level, then continue heating.
- Composed of first distilling, second distilling water and heating two control system.
- Immerge electric heating tube made from copper, heat efficiency is better.



#### Parameters

Model	Electric Supply	Power	Water discharge	Product size (Wooden)	Gross Weight (kg)	Net Weight (kg)
DZ5C	380V/50HZ	10kW	≥5L/h	67 X 43 X 86	32	16
DZ10C	380V/50HZ	16kW	≥10L/h	74 X 43 X 96	38	18
DZ20C	380V/50HZ	32kW	≥20L/h	84 X 52 X 118	56	34