

Automatic Stainless Steel Distiller

DZ5Z / DZ10Z / DZ20Z

Performance and characteristic:

- All of products are made from stainless steel.
- Condenser is made of seamless stainless-steel pipes and so with very good cooling effects & great capacity of water produce.
- The cooled water returning to the vaporizing drum can save energy and water.
- Drain valve can push the water in container, prevent bilge.
- When lack water, will auto stop heating, when the return water reach the normal level, then continue heating.
- Immerge electric heating tube made from copper, heat efficiency is better.



Parameters

Model	Electric Supply	Power	Water discharge	Product size (Paper)	Gross Weight (kg)	Net Weight (kg)
DZ5Z	220V/50HZ	4.5kW	≥5L/h	35X34X78	7	5
DZ10Z	380V/50HZ	7.5kW	≥10L/h	42X38X90	10	8
DZ20Z	380V/50HZ	15kW	≥20L/h	51X36X80	13.5	11.5