

## Common Stainless Steel Distiller

### DZ5 / DZ10 / DZ20

#### Performance and characteristic:

- All of products are made from stainless steel.
- Condenser is made of seamless stainless-steel pipes and so with very good cooling effects & great capacity of water produce.
- The cooled water returning to the vaporizing drum can save energy and water.
- Drain valve can push the water in container, prevent bilge.
- Immerge electric heating tube made from copper, heat efficiency is better.



DZ5

DZ10

DZ20

#### Parameters

Model	Electric Supply	Power	Water discharge	Product size (Paper)	Gross Weight (kg)	Net Weight (kg)
DZ5	220V/50HZ	4.5kW	≥5L/h	35 X 26 X 78	7	5
DZ10	380V/50HZ	7.5kW	≥10L/h	42 X 30 X 90	10	8
DZ20	380V/50HZ	15kW	≥20L/h	51 X 36 X 80	13.5	11.5