

## DW-IID Ultrasonic Homogenizer

As a new type of ultrasonic homogenizer, it enjoys complete functions, novel appearance and reliable performance. It has a large-screen display and concentrated control by a central computer. The ultrasonic time and power can be set accordingly. Besides, such functions as sample temperature display, actual frequency display, computerized tracking and automatic trouble warning can all be displayed in the large LCD screen.



### Parameters

MODEL No	DW-IID
Frequency	20-25KHz
Power	20-1000W Adjustable
Display	7" TFT Touch display
Temperature alarm range	0-99.9C
Timer	1-999 minutes
Accessories (Φ)	Φ6
Working method	Interval /continue
Alarm items	Time, temperature, overload
Optional Tip (Φ)	Φ2、3、10、15