

## JY92-IIDN Ultrasonic Homogenizer



### Parameters

MODEL No	JY92-IIDN
Frequency	20-25KHz
Power	900W (1%- 100%)
Display	7" TFT Touch display
Temperature alarm range	0-99.9C
Timer	0-999 minutes
Accessories (Φ)	Φ 6mm
Working method	Interval /continue
Alarm items	Time, temperature, overload
Optional Tip (Φ)	Φ2mm,3mm,10mm,12mm,15mm
Crushing Capacity	100μl -600ml
Duty Ratio (%)	0.1-99.9%
Supply	220/110Y 50Hz/60Hz
Host size& wight	430mm*245mm*300mm &9.6Kg
Sound chamber size &weight	345mm*345mm*535mm &12kg