

DW-SXT-06 Soxhlet Extractor

DW-SXT-06 Soxhlet Extractor is for labs to make Soxhlet extraction, designed in accordance with the classical Soxhlet extraction method. This apparatus is reasonably designed, provides stable performance, high precision, time-saving operation, is an ideal Soxhlet extraction apparatus for foodstuff, oil&fat, feedstuff and other industries.

Features

- A. Take water bath to heat, providing more even warming;
- B. Solvent for extraction can be recovered and used repeatedly;
- C. Extraction duration can be set up at random, alarm when time is out.

Technical Parameters

- 1.Measurement range: for extract samples like cereals, feedstuff, oily seed, soil, etc.
- 2.Sample number:6/batch
- 3.Capacity of extraction flask:150ml/per
- 4.Sample amount:0-5g/per
- 5.Temp. range:ambient+5°C ~ 100°C
- 6.Power:220V±10V, 50Hz
- 7.Heating power:800W
- 8.Dimensions(mm):550×320×626
- 9.Weight:23kg

