

TU4000E Turbidity Meter

Features

1. 7" colorful touchscreen.
2. Tungsten lamp light source met USEPA 180.1 requirement.
3. 3 photodiodes to ensure high accuracy.
4. 1-8 points calibration using Formazin standards or zero water.
5. Autoranging from 0-4000 NTU.
6. Data Storage capacity of 1000 sets (GLP-compliant).
7. Data stored can be transferred into the computer by USB or RS-232 communication interface.
8. Reset feature automatically resumes all settings back to factory default options.

Includes

- 5 pcs turbidity bottles



Specifications

Model		TU4000E
Measuring Unit	Range	(0.000~20.00) NTU, (20.01~200.0) NTU, (200.1~2000) NTU, (2001~4000) NTU
	Accuracy	±6%
	Repeatability	≤0.5%
	Fluctuation	±0.5% FS/30min
General	Power Supply	AC Adapter, 200-240 VAC input
	Dimensions	358×323×160 mm
	Weight	6.5 kg
	Dimensions (Carton)	490×410×220 mm
	Weight (Gross)	7.5 kg