

COD-100R Reactor

Features

1. 21 samples at a time.
2. Settable heating temperature between 100-165 °C.
3. Settable heating time between 0-120 min.
4. Transparent plastic protection shield is equipped to protect the operator.

Includes

- Test-tube rack
- 21 pcs digestion tubes



Specifications

Model		COD-100R
Measuring Unit	Temperature	100-165 °C
	Time	0-120 min
General	Power Supply	AC Adapter, 200-240 V AC input
	Dimensions	330×228×132mm
	Weight	4 kg
	Dimensions (Carton)	575X435X285 mm
	Weight (Gross)	5.5 kg