

DW-ZBJ01

Intelligent Disintegration Apparatus

Product description

DW-ZBJ01 intelligent disintegration apparatus is a new type of apparatus designed in accordance with the requirements of the Chinese Pharmacopoeia. The detection of the disintegration time limit of 18 samples in 3 cups can be achieved.

Features

1. 3 Basket Disintegration Time Limit Detection System
2. Adjustable test time
3. Real-time display of time and temperature
4. Overheating protection
5. Integral water bath box, no dead angle of water bath circulation, easy to clean
6. Overall appearance design, high-grade and concise
7. Independent lighting system
8. High sensitivity capacitive touch screen
9. Authority management
10. Audit trail.



Parameters

| | |
|--|--|
| Model | DW-ZBJ01 |
| Temperature range | Room temperature~45°C |
| Temperature Resolution | 0.01°C |
| Number of hanging baskets | 3 cups |
| Lifting frequency | 30-32 times/min |
| Lifting distance | 55mm±1mm |
| The distance between the bottom of the screen and the bottom of the cup | 25mm±2mm |
| Mesh aperture | 2mm, optional 1mm, 0.71mm, 0.42mm |
| Block density | 1.18—1.20g/cm³ |
| Number of users | 7 users |
| User level | III |
| Record storage | 96*366 days |