

DW-E/F Weighing Scale

Features

- Stable weighing system
- LED display
- Full range tare & Counting, units conversion functions
- Set mini weighing by users, and choose the suitable interface
- Stainless steel large scale with sprayed carbon steel base
- Optional RS232 can connect with PC and other instruments



Parameters

Model	XY60E	XY100E	XY150E	XY200E	XY300E	XY500E	XY60F	XY100F	XY150F
[kg]Max capacity	61	110	160	210	310	510	61	110	160
[g]Resolution	10					50	1		
[g]Repeatability	±30					±150	±3		
[g]Linearity	±50					±250	±5		
Resp. Time	3sec								
Temp Range	+5...+35°C								
Sensitivity Drift	±3ppm/°C								
Pan size	400*500mm/450*600mm					600*800mm	400*500mm		
Line Voltage	220V/60Hz								