

DWT-B Series Pulsating Vacuum Table Top Steam Sterilizer

Features

1. Quality stainless steel inner tank, sealed cover
2. Safety interlocking device
3. Automatic water
4. With drying function
5. With BD test, HELIX test, VACUUM test.
6. LCD display and touch buttons.
7. With standard test interface
8. European CLASS B standard
9. With automatic protection function: over temperature protection, over pressure protection, low water level protection, anti dry burning, and safety door lock.
10. Three times per-vacuum
11. Buzzer reminder after sterilization, automatic stop.
12. Automatic discharge of cold air, automatic steam exhaust after sterilization.
13. Built-in steam generator to provide saturated steam quickly
14. With printer,USB function and U disk
15. Microcomputer control

DWT-8B/12B



DWT-18B/23B

Parameters

Mode	DWT-8B	DWT-12B	DWT-18B	DWT-23B
Volume	8L	12L	18L	23L
Power	2kw	2kw	2kw	2kw
Voltage	110V/220V, 60Hz/50Hz			
Design pressure	0.25MPa			
Design temperature	139°C			
Rated working pressure	0.22MPa			
Rated working temperature	134°C			
Chamber volume	Φ170×320	Φ200×355	Φ250×370	Φ250×452
Outer dimension	445×580×400	445×580×400	480×640×450	480×690×450
Packing size(LxWxH)	650×500×450	650×500×450	720×560×510	760×560×510
Gross/net weight	45/42kg	48/44kg	57/48kg	60/53kg