

DWT-BA Series Pulsating Vacuum Table Top Steam Sterilizer

Features

1. Automatic door switch
2. Automatic water
3. Quality stainless steel inner tank, sealed cover
4. Microcomputer control
5. Safety interlocking device
6. With drying function
7. With BD test, HELIX test, Vacuum test.
8. LED display and inductive operation
9. Equipped with standard test interface
10. European CLASS B standard.
11. With automatic protection function: over temperature protection, over pressure protection, low water level protection, anti dry burning, safety door lock.
12. Three times per-vacuum
13. Buzzer reminder after sterilization, automatic stop.
14. Automatic discharge of cold air, automatic steam exhaust after sterilization.
15. Built-in steam generator to provide saturated steam quickly
16. With built-in printer, USB function and U disk



Parameters

Mode	DWT-18BA	DWT-23BA	DWT-45BA
Volume	18L	23L	45L
Power	2kw	2kw	2.5kw
Voltage	110V/220V, 60Hz/50Hz		
Design pressure	0.25MPa		
Design temperature	139°C		
Rated working pressure	0.22MPa		
Rated working temperature	134°C		
Chamber volume	Ø250×370	Ø250×470	Ø316×615
Outer dimension	620×505×450	505×710×450	657×870×575
Packing size(LxWxH)	670×560×510	760×560×510	980×750×630
Gross/net weight	62/53kg	65/58kg	142/105kg