

DW-280B Electric Coal and Gas Heating Portable Autoclave

Features

1. High Quality SUS304 stainless steel material
2. Electric, coal and gas heating
3. Silicone seal
4. Dual scale indicate pressure gauge
5. Double safety valve to ensure safe
7. Easy operation, safe and reliable



Specification

Model	DW-280B	DW-280B	DW-280B
Volume	18L	24L	30L
Chamber size	Φ300*330mm	Φ300*460mm	Φ300*520mm
Working temperature	115-129°C	115-129°C	115-129°C
Working pressure	0.145-0.165Mpa	0.145-0.165Mpa	0.145-0.165Mpa
Packing size	40*40*42cm	40*40*53cm	40*40*64cm
G.W/N.W	16kg/15kg	17kg/16kg	18kg/117kg
Power	AC220V/50Hz/2.0KW	AC220V/50Hz/2.0KW	AC220V/50Hz/2.0KW