

DW-280D Portable Steam Autoclave

Features

1. High Quality SUS304 stainless steel material
2. Electric heating system
3. Silicone seal
4. Dual scale indicate pressure gauge
5. Double safety valve to ensure safe
6. High pressure faucet drainage water
7. Over-temperature protection
8. Over-pressure protection
9. Water shortage protection
10. Easy operate, safe and reliable



Specification

Model	DW-280D	DW-280D	DW-280D
Volume	18L	24L	30L
Chamber size	Φ300*330mm	Φ300*460mm	Φ300*520mm
Working temperature	115-129°C	115-129°C	115-129°C
Working pressure	0.145-0.165Mpa	0.145-0.165Mpa	0.145-0.165Mpa
Time scope	0-80min	0-80min	0-80min
Packing size	40*40*53cm	40*40*64cm	40*40*70cm
G.W/N.W	18kg/16kg	19kg/17kg	20kg/18kg
Power	AC220V/50Hz/2.0KW	AC220V/50Hz/2.0KW	AC220V/50Hz/2.0KW