

DT Series Ultrasonic Cleaner

Product Introduction

Ultrasonic cleaning machine is to use the cavitation effect of ultrasonic waves in the cleaning liquid to accelerate the dissolution and stripping of dirt, grease and other dirt in the inner hole and blind hole of the vessel being cleaned. DRAWELL ultrasonic cleaning machine adopts imported key equipment to ensure that the instrument has the advantages of fast cleaning speed, good effect, no damage to the surface and reduced labor intensity, etc.



Parameters

Model	Inner tank (mm)L/W/H	Volume (L)	Frequency (KHz)	Ultrasonic power (W)	Heating power (W)	Temperature adjustable (°C)	Time adjustable (min)	Drainage	Rack
DW-100DT	300*150*100	4.5	40	100	200	RT-80	1-99	/	√
DW-120DT	240*140*100	3	40	120	200	RT-80	1-99	/	√
DW-120DT	240*140*150	5	40	120	200	RT-80	1-99	/	√
DW-3200DT	300*155*150	6	40	180	400	RT-80	1-99	/	√
DW-5200DT	300*240*150	10	40	240	800	RT-80	1-99	√	√
DW-5200DT	300*240*150	10	40	360	800	RT-80	1-99	√	√
DW-4200DT	300*240*200	14.4	40	480	800	RT-80	1-99	√	√
DW25-12DT	500*300*150	22.5	40	600	1000	RT-80	1-99	√	√
DW25-12DT	500*300*150	22.5	40	720	1000	RT-80	1-99	√	√
DW-800DT	500*300*200	30	40	840	1000	RT-80	1-99	√	√
DW-1000DT	600*300*300	54	28	1000	3000	RT-80	1-99	√	√
DW-1200DT	600*400*300	72	28	1200	4000	RT-80	1-99	√	√
DW-1500DT	600*500*300	90	28	1500	4000	RT-80	1-99	√	/
DW-1800DT	600*600*300	108	28	1800	4000	RT-80	1-99	√	/
DW-2000DT	650*600*300	117	28	2000	6000	RT-80	1-99	√	/
DW-3600DT	750*600*300	135	28	3600	8000	RT-80	1-99	√	/
DW-6000DT	1200*650*350	273	28	6000	10000	RT-80	1-99	√	/