

## DL-7MC

# Inelligent Type Large Capacity Refrigerated Centrifuge

### Features:

1. With the self-diagnostic and protection system for imbalance, over-speed, over-temperature. Greatly ensured the safe of users and machine.
2. Electric control lid locks, can control both close and open door lid, with backup power to avoid power off unexpected, Backup power automatically recharges when machine powered on.
3. With large 7 Inches color panel, can operate with gloves.
4. With wheels easy to short distance move, and wheels can be fixed.
5. Automatic Rotor Identification to avoid over-high speed running of rotor.
6. Data can be stored (Option function).



### Parameters

Model	DL-7MC	Max Speed	8000r/min
Max RCF	11620xg	Max Capacity	6×1000ml
Speed Accuracy	±10r/min	Timer Range	0-99h 59 min
Temperature Range	-20~+40°C	control&drive	microprocessor control/direct drive
Acc/dec Rate	0-9(Grade)	Program Store	99
Motor Power	4KW	Refrigerator Power	2.2KW
Noise	≤60dB(A)	Power Supply	AC220V/380V 50Hz 35A
Weight(Without Rotor)	280kg	Dimensions	730×840×100mm(L×W×H)
Refrigeration System	Big Brand High-performance Compressor CFC-free Refrigeration System		
Rotor	No.1 AR 4*1000ml (8000rpm 11620xg) No.2 AR 6*1000ml (7000rpm 11500xg)		