

## Automatic pH Analyzer

### Features


1. Touch screen, back light adjustable, parameter setting.
2. Real-Time clock.
3. Higher efficiency: up to 40 sample solution in a single batch.
4. GLP data management: samples ID classification.
5. Reliable test: measurements traceable to raw data.
6. Compact design, highly optimize the test working room.
7. Automated electrode cleaning and blowing dry.
8. Durable and reliable aluminum alloy tray, easy to clean, corrosive-proof.

Conforms to standards of: GB/T 7573-2009

Application: Textiles, Leather samples/extraction solution



Working station pan view-Aluminum alloy



### Specifications

<b>pH range</b>	<b>-2-19.99pH (± 0.02pH ATC)</b>
<b>Calibration</b>	<b>1- 3 points, automatically</b>
<b>Measurement</b>	<b>automatically</b>
<b>Purified water refilling</b>	<b>automatically</b>
<b>Wastewater collection</b>	<b>automatically</b>
<b>GLP management</b>	<b>sample IDs setting</b>
<b>Temperature range</b>	<b>0 – 100.0 °C(±0.5°C)</b>
<b>Data storage</b>	<b>2000 readings</b>
<b>Data export</b>	<b>PC/Flash Disk</b>
<b>Power supply</b>	<b>220V, 50Hz, 1000W</b>
<b>Meter size</b>	<b>78.0×42.5×45.5cm 26kg</b>