

Floortype High Speed Refrigerated centrifuge GL-10MC

Features:

Superior performance

Microcomputer controlled, AC frequency variable motor drive, without power pollution, and longer lifetime. Lower noise and less vibration, equipped with multiple types of rotors for users to choose. Efficiency has been greatly improved.

Easy operation

Equipped with Large LCD with touch function display, programmable application, which can be operated freely with variety of rotors for selection.

Good effect

With function of freely setting the acceleration time, deceleration time, and automatically calculating RCF, which can achieve good centrifugal effect.

Safe and reliable performance

With protection from over speed, over temperature, imbalance and protection of automatic electronic or mechanical lock, to ensure the safety of operators and machines.



Model	GL-10MC	
Max Speed/RCF	10000r/min 18300xg	
Speed Accuracy	±20r/min	
Temp Range	-20°C~+40°C±1°C	
Refrigeration system	High-performance Compressor CFC-free Refrigeration System	
Timer Range	0~23h 59min	
Motor	Microprocessor control, Frequency conversion motor	
Noise	≤60dB(A)	
Power Supply	AC220V 50-60Hz 30A	
Dimension	$710\times820\times1200$ mm(L×W×H)	
Net Weight	260kg	

Rotor No.	Max Capacity	Max Speed	Max RCF
No.1 AR	6x300ml	10000r/min	18300xg
No.2 AR	6x500ml	8000r/min	11620xg
No.3 SR	4x1000ml	4000r/min	4050xg