

SHZ-95B Water Circulation Vacuum Pump

Features:

- ◆ 5 exhaustors can be used individually or in parallel, with a vacuum gauge.
- ◆ Pump body is made of stainless steel; motor made by Sino-foreign joint venture; pump spindle made of stainless steel; pump head design complying with fluid mechanics; shell is made of engineering plastics or stainless steel.
- The exhaust speed is as fast as 14L/MIN.
- Anti-corrosive, without pollution, low noise, easy to move.
- Optional: vacuum regulating valve.
- This vacuum pump can absorb the harmful gas, which is environmental friendly and ensures the physical and psychological health of experimenters.
- Water and energy saving instrument.



Parameter

MODEL #	SHZ-95B
Power (W)	370
Voltage (V/Hz)	220/50
Flow Rate (L/min)	≥80
Lift (m)	≥12
Safety Protection	Check valve
Shell material	Cold rolled sheet with anti - static spray
Max. Vacuum Degrre (MPa)	0.098
Unit Extraction (L/min)	10
Exhaustor (A)	5
Water tank volume (L)	≥50
Water tank material	PVC
Product size (mm)	450L×350W×820H