

SZ Series

Automatic Pure Water Distiller

Features:

- ◆ 1. The SZ series automatic pure water distiller is heated with quartz heating tube, which is energy-saving, pure, durable and of high radiation rate.
- ◆ 2. The adoption of stainless steel in supporting shelf and clamp; special engineering plastics made of control panel and mian distiller body made of quality glass.
- ◆ 3. The pure water distillated does not contact with any metal products.
- 4. The SZ series automatic water distiller is ideal for the distillation of pure water.
- 5. With different requirement of the pure water, SZ series automatic water distiller can be classified into the following 3 types:

Single stage automatic pure water distiller	Double stage automatic pure water distiller	Triple stage automatic pure water distiller
Pure water	Purer water	Purest water

- 6. All SZ series automatic pure water distiller is equipped with reply protection, which is safe and reliable to use.
- ◆ 7. The distiller switch can be controlled with temperature controller and spring gauge controller, which helps avoid glass breakage caused by dry heating of loop heat pipe.
- 8. Power supply will be automatically connected or cut off in high water line or low water line. So the pure water distiller is fully automatic.

Parameter

MODEL #	SZ-96	SZ-93	SZ-97
Flow rate	1800ml/h	1600ml/h	1500ml/h
Input power	1.5Kw	3Kw	4.5Kw
Voltage	AC220v/50Hz	AC220v/50Hz	AC220v/50Hz
Pack size LWH(cm)	Pack size LWH(cm)	Pack size LWH(cm)	Condenser 57*48*39cm,5.8KG Host 61*45*44cm,10KG
Weight	5.4Kg	10Kg	15.8Kg
Size (cm)	77×20×68cm	77×36×78cm	77×49×89
Picture			