

LCD Digital Rotator v.s Economical Rotator

Model: MX-RD-Pro / MX-RD-E

Rotators give gentle but effective mixing of biological samples held in 1.5mL to 50mL tubes, and are used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics, etc.

MX-RD-Pro Features:

- The rotator comes with a disk type design.
- LCD displays time and speed.
- Adjustable speed range from 10 to 70rpm.
- Adjustable mixing angle between 0 and 90 o
- Used with a wide variety of accessories
- Continuous and timed quiet operation

MX-RD-E Features:

- Provides gentle and effective mixing.
- Adjustable speed range of 0 to 80rpm
- Used with a wide variety of accessories
- Continuous quiet operation





Parameters

Model	MX-RD-Pro	MX-RD-E	
Motor type	DC motor	DC motor	
Tilt angle	0-90°		
Speed range	10-70rpm	0-80rpm	
Speed display	LCD	Scale	
Timer display	LCD	Scale	
Timer	1-1199min		
Operation mode	Continuous / timer	Continuous	
Voltage	100-240V, 50/60Hz		
Power	40W	20W	
Dimension [D×W×H]	280×210×300mm	150X530X190mm	
Permissible ambient temperature and humidity	5-40oC, 80%RH		
Protection class	IP21	IP21	



Accessory

Cat. No.	Picture	Descriptions	
18900160		Disk accessory, for 1.5ml x 60 centrifuge tubes holder, use with MX-RD-Pro	
18900161		Disk accessory, 15ml x 16 centrifuge tubes holder, use with MX-RD-Pro	
18900162		Disk accessory, 50ml x 8 centrifuge tubes holder, use with MX-RD-Pro	
18900140		Disk support rods, 4 pcs, use with disk accessories	
18900327		Disk and Rotisserie accessory, for 1.5mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E	
18900328		Disk and Rotisserie accessory, for 15mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E	
18900329		Disk and Rotisserie accessory, for 50mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E	
18900196		1.5mL metaltube holder,use with MX-RL/RD-Pro	
18900197		15mL metal tube holder price, use with MX-RL/RD-Pro	
18900198		50mL metal tube holder price, use with MX-RL/RD-Pro	