

# LCD Digital Rotator v.s Classic Rotator

## Model: MX-RL-Pro / MX-RL-E


The rotators give gentle but effective mixing of biological samples to be held in 1.5ml to 50ml micro tubes, used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics etc.


### MX-RL-Pro Features:

- ◆ Electronic adjustment of speed and time
- LCD displays speed and time
- Speed range from 10 to 70rpm
- Provides gentle but effective rolling
- A wide variety of accessories are available

### MX-RL-E Features:

- Provides gentle but effective rolling.
- Adjustable speed range from 0 to 80 rpm.
- ♦ A wide variety of accessories are available.
- Continuous and timed quiet operation.





#### Parameters

Model	MX-RL-Pro	MX-RL-E
Motor type	DC motor	DC motor
Speed range	10-70rpm	0-80rpm
Speed display	LCD	Scale
Timer display	LCD	Scale
Timer	1-1199min	
Operation mode	Continuous/timer	Continuous
Voltage	100-240V, 50/60Hz	100-240V, 50/60Hz
Power	<b>40W</b>	20W
Dimension[D×W×H]	220×510×260mm	300×220×310mm
Permissible ambient temperature and humidity	5-40oC,80%RH	5-40oC, 80%RH
Protection class	IP21	IP21



#### Accessory

Cat. No.	Picture	Descriptions	
18900140		Disk support rods, 4 pcs, use with disk accessories	
18900142		Rotisserie accessory, for 1.5ml x 48 centrifuge tubes held horizontally, use with MX-RL-Pro	
18900143		Rotisserie accessory, for 15ml x 24 centrifuge tubes held horizontally, use with MX-RL-Pro	
18900144		Rotisserie accessory, for 50ml x 24 centrifuge tubes held horizontally, use with MX-RL-Pro	
18900145		Rotisserie accessory, for 1.5ml x 32 centrifuge tubes held vertically, use with MX-RL-Pro	
18900146		Rotisserie accessory, for 15ml x 16 centrifuge tubes held vertically, use with MX-RL-Pro	
18900147		Rotisserie accessory, for 50ml x 16 centrifuge tubes held vertically, use with MX-RL-Pro	
18900327		Disk and Rotisserie accessory, for 1.5mlx 12 centrifuge tubes,	
	use with MX-RL-E and MX-RD-E PC		
18900328		Disk and Rotisserie accessory, for 15mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E PC	
	3900329	Disk and Rotisserie accessory, for 50mlx 12 centrifuge tubes,	
18900329		use with MX-RL-E and MX-RD-E PC	
18900196、18900197、18900198		1.5mL 15mL 50mLmetaltube holder, use with MX-RL/RD-Pro	