



DW-C High Pressure Vertical Autoclave

◆ Features

- SUS 304 stainless steel sterilizing chamber
- Accurate dual scale indicate pressure gauge
- Silicone seal ensuring no steam leakage
- Spin buttons to set sterilization temperature and time
- Blue, orange and red light showing work status
- Automatic release when overpressure
- Water shortage and over heating protection
- Easy operate, safe and reliable



◆ Specification

Model	DW-C35L	DW-C50L	DW-C75L
Volume	35L	50L	75L
Chamber size	Φ300*485mm	Φ320*545mm	Φ388*480mm
Working temperature	115-129℃	115-129℃	115-129℃
Working pressure	0.145-0.165Mpa	0.145-0.165Mpa	0.145-0.165Mpa
Time scope	0-80min	0-80min	0-80min
Packing size	53*52*92cm	56*56*105cm	62*60*109cm
G.W/N.W	49kg/34kg	60kg/40kg	80kg/57kg
Power	AC220V/50Hz/2.5KW	AC220V/50Hz/2.5KW	AC220V/50Hz/3.5KW

