



DW-C Digital Vertical Autoclave

◆ Features

- SUS 304 stainless steel sterilizing chamber
- Accurate dual scale indicate pressure gauge
- Silicone seal ensuring no steam leakage
- LED digital display temperature and time of sterilization
- Temperature and time of sterilization can be set
- Automatic release when overpressure
- Water shortage and over heating protection



◆ Specification

Model	DW-C Digital	DW-C Digital	DW-C Digital
Volume	35L	50L	75L
Chamber size	Φ300*485mm	Φ320*545mm	Φ388*480mm
Working temperature	100-134℃	100-134℃	100-134℃
Working pressure	0.22Mpa	0.22Mpa	0.22Mpa
Time scope	0-999min	0-999min	0-999min
Packing size	53*52*92cm	56*56*105cm	62*60*109cm
G.W/N.W	49kg/34kg	60kg/40kg	80kg/57kg
Power	AC220V/50Hz/2.5KW	AC220V/50Hz/2.5KW	AC220V/50Hz/3.5KW

