



Hot(micro/poarization) Series

Improve the efficiency of innovative points

Closed furnace, guarantee the stability of temperature;
Accurately the heating rate of 0.1 to 20.0 and infinitely adjustable;
TFT true color touch-screen;
Can be customized as many as 100 kinds of plan;
Measurable block ,paste material;
Can connect printers, print data in real-time.

◆ The Main Technical Indicators

Model	JH350	JH370
Temperature range	RT-300°	RT-350°
Temperature resolution	0.1℃	0.1℃
Temperature gradients	(0.1-20.0) ℃/Min	(0.1-20.0) ℃/Min
Temperature accuracy	± 0.3℃	± 0.3℃
Slide/cover glass specification	76*19*1/15*15*0.5mm	76*19*1/15*15*0.5mm
Field of vision	2mm	2mm
X-y table	Most mobile 13 mm	Most mobile 13 mm
Working distance of lens	At least 7 mm	At least 7 mm
Result storage	200	400
Display	TFT Touch screen	TFT Touch screen
Interface	USB RS232 SD card Enetnet	USB RS232 SD card Enetnet
Power	110-220V 50/60HZ 120W	110-220V 50/60HZ 120W

