

Fixed Top LCD Digital Orbital Shaker

Model: SK-O330-M

Orbital Shaker is ideal for applications like bacterial suspension, solubility studies, growth of bacteria and yeast, extraction procedures, washing procedures, diagnostics tests, hybridization, general mixing, staining & destaining, mixing, and extracting.

Features

- Max. load capacity 3kg;
- Independent LCD display for speed and time;
- Speed and time value can be altered during running;
- Maintenance-free brushless DC motor;
- Foot screws height can be adjusted to make sure stable running;
- Continuous or timing mode optional.



Specifications

Model#	SK-O330-M
Orbital diameter	10mm
Max. load capacity (with platform)	3kg
Motor type	Brushless DC motor
Motor input power	20W
Motor output power	15W
Power	25W
Time settings range	1min-19h59min
Speed range	70~400rpm
Speed display	LCD
Timer display	LCD
Voltage	100-240V,50/60Hz
Protection class	IP21
Dimension [W x D x H]	420x370x100mm
Weight	9.8kg
Permissible ambient temperature and humidity	5-40°C 80%RH