

Tower Electric Heating Stainless Steel Distiller

TZ50/TZ100/TZ200/TZ400/TZ600

Performance and characteristic:

- All of products are made from stainless steel.
- Use high temperature steam to heat, save the energy.
- Use the steam condensation water.
- Seamless stainless-steel pipes, heat efficiency is better.
- Condenser is made of seamless stainless-steel pipes and so with very good cooling effects & great capacity of water produce.
- Percolation, pushing ammonia, dividing water and steam when distilling, in order to control the quality.



Parameters

Model	Water discharge	Pressure of steam supply	Pressure of cooling water	Steam consumption	Ratio of water consumption	Product size (Wooden)	Gross Weight (kg)	Net Weight (kg)
TZ50	≥50L/h	0.4MPa	0.2MPa	≈1:1	≤1:14	75 X65 X145	88	55
TZ100	≥100L/h					77 X 67X 165	100	65
TZ200	≥200L/h					102 X 82X140	150	95
TZ400	≥400L/h					113 X103 X165	200	130
TZ600	≥600L/h							