

## JY92-IIN Ultrasonic Homogenizer



### Parameters

| MODEL No                   | JY92-IIN                    |
|----------------------------|-----------------------------|
| Frequency                  | 20-25KHz auto tracking      |
| Power                      | 650W (1%- 100%)             |
| Display                    | 7" TFT Touch display        |
| Temperature alarm range    | 0-99.9C                     |
| Timer                      | 0-999 minutes               |
| Accessories (Φ)            | Φ6mm                        |
| Working method             | Interval /continue          |
| Alarm items                | Time, temperature, overload |
| Optional Tip (Φ)           | Φ2mm, 3mm, 10mm             |
| Crushing Capacity          | 100μl -500ml                |
| Duty Ratio (%)             | 0.1-99%                     |
| Supply                     | 220/110Y 50Hz/60Hz          |
| Host size& wight           | 430mm*245mm*300mm &9.5Kg    |
| Sound chamber size &weight | 345mm*345mm*535mm &12kg     |