

JY98-IIIDN Ultrasonic Homogenizer



Parameters

MODEL No	JY98-IIIDN
Frequency	19.5-20.5KHz
Power	1200W (1%- 100%)
Display	7" TFT Touch display
Temperature alarm range	0-99.9C
Accessories (Φ)	Φ 20mm
Working method	Interval /continue
Alarm items	Time, temperature, overload
Optional Tip (Φ)	Φ15mm, 25mm
Crushing Capacity	50- 1000ml
Duty Ratio (%)	0.1-99.9%
Supply	220/110V 50Hz/60Hz
Host size& wight	430mm*245mm*300mm &11.6Kg
Sound chamber size &weight	345mm*345mm*535mm &12kg