

DW- KFT-40C Karl Fischer Titrator

Features

- LCD display screen.
- Support KF Titration Mode and Titer Detection Mode.
- Features of auto-filling, auto-purging and auto-mixing of the reagents ensure safe handling of Karl Fischer chemicals.
- Selectable units including mg, mg/L, %, ppm, etc.
- Settable parameters, including measurement unit, polarization current, stirring rate, titration rate, stop volume, endpoint potential, stop criterion and etc.
- Store up to 200 titration data sets (GTP-compliant).
- Data can be easily transferred to printer via RS-232 commnication interface.
- Reset feature automatically resumes all settings back to factory default options.

Includes

1. Measuring electrode
2. 10ml burette
3. Volumetric titration cell



Specifications

Model		DW-KFT-40C
Volumetric Titration	Range	0.1mg~250mg
	Resolution	0.1mg
	Working Current	1μA±0.2μA; 50μA±10μA
	Repeatability	≤0.5%
General	Power Supply	AC Adapter, 100-240 V AC input
	Dimensions	340x400x400 mm
	Weight	10kg
	Dimensions (Carton)	540x490x480 mm
	Weight (Gross)	13 kg