

DW-Ti-20 Potential Titrator

Features

- LCD display screen.
- Support the following titration modes: SET (Preset Endpoint Titration) and MAT (Manual titration).
- Stirrer is included.
- Glass burette can be used for volume measurement (prepared by user).

Includes

1. JB-10 stirrer
2. 217-01 reference
3. 213-01 platinum
4. 216-01 silver ele
5. E-201F pH comp
6. NIST pH standar



Specifications

Model		DW-Ti-20
Measuring Unit	Range	(0~±1400) mV, (0~14.00)pH
	Resolutions	1mV, 0.01pH
	Accuracy	pH: ±0.03pH, mV: ±0.35%FS
	Fluctuation	±0.01pH/3h
General	Power Supply	AC Adapter, 200-240 VAC input
	Dimensions	300×235×100 mm
	Weight	5 kg
	Dimensions (Carton)	370×305×330 mm
	Weight (Gross)	6 kg