

EC110B Conductivity Meter

Features

- LCD display screen, 3.5 inches.
- Reset feature automatically resumes all settings back to factory default options.
- Auto-power off feature effectively extends the battery service life.
- IP65 waterproof. The portable meter is suitable for fields measurements and out door measurements.
- Cell constant is settable
- Manual temperature compensation ensures accurate results.

Includes

- ★ DJS-1VC EC electrode
- ★ 1413μS/cm EC standard solution, 50ml, 1 vial
- ★ Electrode holder
- ★ Silicone rubber case
- ★ Wristbands



Specifications

Model		EC110B
Parameters		EC
Conductivity	Range	0.00 μS/cm to 100 mS/cm
	Resolution	0.01 μS/cm minimum, various with range selection
	Accuracy	±1.0% FS
	Reference Temperature	25 °C
Measurement	Reading Mode	Continuous
	Reading Prompts	Reading, Stable
	Temp. Compensation	MTC
Inputs	Conductivity	5-pin aviation connector
Display options	Backlight	Yes
	Auto Shutdown	300, 600, 1200, 1800, 3600 sec., off
	IP Rating	IP65
General	Power	Rechargeable Lithium battery, AC Adapter, 100-240Y AC input, DC5V output
	Dimensions	80×225×35 mm
	Weight	400g (0.88 lb)
	Dimensions (Carton)	300×290×140 mm
	Weight (Gross)	1 kg