

DW-MW Series、DW-MW-A Halogen Moisture Analyzer

Features

- ◆ Halogen lamp/HBM sensor
- ◆ Aluminum shell/Stainless steel chamber
- ◆ LCD(White back light with black font)
- ◆ Temperature/Time set
- ◆ The percentage of moisture content/dry residual
- ◆ Store historical set
- ◆ Option:Aluminum box/Printer/Grinder/USB



DW-MW

- ◆ Halogen lamp/HBM sensor
- ◆ Aluminum shell/Stainless steel chamber
- ◆ High definition LCD
- ◆ Highlight halogen lamp
- ◆ Temperature/Time set
- ◆ The percentage of moisture content/dry residual
- ◆ Store historical set
- ◆ Option:Aluminum box/Printer/Grinder



DW-MW-A



Computer terminal



Printer

Specifications

Model	DW-110MW	DW-105MW	DW-102MW	DW-100MW	DW-100MW-A
Capacity	110g				
Readability	0.01g	0.005g	0.002g	0.001g	
Operating Temp	5°C-35°C				
Display	LCD(White back light with black font)				High definition LCD
Calibration	External Calibration				
Interface	RS232 (standard) / USB (optional)				
Pan Size	Φ90mm				
Packing Size	490x350x360mm				
Gross Weight	7.5kg				
Heating Source	Halogen Lamp-1				Halogen Lamp-2
Temp.Sensor	PT100				
Temp.Range	40°C-199°C				
Temp,Step	1°C				
Moisture Range	0.09%-100.0%	0.009%-100.00%			
Moisture Readability	0.2%	0.10%	0.04%	0.01%	
Dry Rang	100.00%-0.00%				
Dry Readability	0.2%	0.10%	0.04%	0.01%	
Temp Setting	40°C-199°C By 1°C Step				
Time Setting	1-99min By 10s				
Number of Storage	15				
Heater	220V±15% 50HZ/110V±15% 60HZ/400W				
Balance	Input: 220V 50HZ/110V 60HZ Output:9V				
Chamber Height	22mm				