

DW-B Hand Wheel&Digital Vertical Autoclave

Features

- 1. Full SUS 304 stainless steel
- 2. PID Microcomputer Control heating, air release, temperature, time and steam realease
- 3. LED digital display sterilization temperature and time
- 4. Eland wheel quick open door
- 5. Double safety interlocking system to ensure safe
- 6. Automatic release when overpressure
- 7. Water shortage and over heating protection
- 8. Tempered glass anti-scald protection
- 9. Sound-light alarm when sterilization finished
- 10. Drying function can be addes if customer need (best more than 75L)



Specification

Model	DW-B hand wheel&digital	DW-B hand wheel&digital	DW-B hand wheel&digital	DW-B hand wheel&digital
Volume	35L	50L	75L	100L
Chamber size	Φ 318*420mm	Φ 318*530mm	Φ 388*480mm	Φ 388*630mm
Working temperature	100-134℃	100-134℃	100-134℃	100-134℃
Max Working pressure	0.22Mpa	0.22Mpa	0.22Mpa	0.22Mpa
Time scope	0-9909min	0-9909min	0-9999min	0-9999min
Packing size	55*54*100cm	57*55*109cm	62*62*121cm	62*62*131cm
G.W/N.W	65kg/44kg	73kg/49kg	91kg/71kg	103kg/74kg
Power	AC220V/50Hz/2.5KW		AC220V/50Hz/3.5KW	AC220V/50Hz/4.0KW