

DW-B Digital Vertical Autoclave

Features

- 1. High quality SUS304 stainless steel material
- 2. Dual scale indicate pressure gauge
- 3. Light indicate working status
- 4. LED digital display screen
- 5. Microcomputer control system
- 6. Temperature, time of sterilization can be set.
- 7. Water lack auto-protection
- 8. Over pressure auto-discharging
- 9. Over current auto-protection
- 10. Over temperature auto-protection
- 11. Drying function can be add(best large than 75L)



Specification

Model	DW-B35LDigital	DW-B50L Digital	DW-B75L Digital	DW-B100L Digital	DW-B120L Digital	DW-B150L Digital
Volume	35L	50L	75L	100L	120L	150L
Chamber size	Φ300*485mm	Φ 320*545mm	Φ 388*480mm	Φ 388*630mm	Φ 460*590mm	Φ 460*700mm
Working temperature	100-134℃	100-134℃	100-134℃	100-134℃	100-134℃	100-134℃
Working pressure	0.22Mpa	0.22Mpa	0.22Mpa	0.22Mpa	0.22Mpa	0.22Mpa
Time scope	0-999min	0-999min	0-999min	0-999min	0-999min	0-999min
Packing size	53*52*92cm	56*56*105cm	62*60*109cm	63*61*118cm	66*67*121cm	67*68*128cm
G.W/N.W	49kg/34kg	60kg/40kg	80kg/57kg	85kg/66kg	97kg/75kg	110kg/80kg
Power	AC220V/50Hz/ 2.5KW	AC220V/50Hz/ 2.5KW	AC220V/50Hz/ 3.5KW	AC220V/50Hz/ 4.0KW	AC220V/50Hz/ 4.0KW	AC220V/50Hz/ 4.0KW