

DW-B Vertical Autoclave

Features

1. High quality SUS304 stainless steel material
2. Dual scale indicate pressure gauge
3. Light indicate working status
4. Temperature, time of sterilization can be set.
5. Water lack auto-protection
6. Over pressure auto-discharging
7. Over current auto-protection
8. Over temperature auto-protection
9. Drying function can be add(best large than 75L)



Specification

Model	DW-B35L	DW-B50L	DW-B75L	DW-B100L	DW-B120L	DW-B150L
Volume	35L	50L	75L	100L	120L	150L
Chamber size	Φ300*485mm	Φ320*545mm	Φ388*480mm	Φ388*630mm	Φ460*590mm	Φ460*700mm
Working temperature	115-129°C	115-129°C	115-129°C	115-129°C	115-129°C	115-129°C
Working pressure	0.145-0.165Mpa	0.145-0.165Mpa	0.145-0.165Mpa	0.145-0.165Mpa	0.145-0.165Mpa	0.145-0.165Mpa
Time scope	0-80min	0-80min	0-80min	0-80min	0-80min	0-80min
Packing size	53*52*92cm	56*56*105cm	62*60*109cm	63*61*118cm	66*67*121cm	67*68*128cm
G.W/N.W	49kg/34kg	60kg/40kg	80kg/57kg	85kg/66kg	97kg/75kg	110kg/80kg
Power	AC220V/ 50Hz/2.5KW	AC220V/ 50Hz/2.5KW	AC220V/ 50Hz/3.5KW	AC220V/ 50Hz/4.0KW	AC220V/ 50Hz/4.0KW	AC220V/ 50Hz/4.0KW