

DW-280 Digital Portable Autoclave

Features

- I. High Quality SUS304 stainless steel material
- 2.. Electric heating system
- 3. Silicone seal
- 4. Dual scale indicate pressure gauge
- 5. Double safety valve to ensure safe
- 6. LED digital display screen
- 7. Microcomputer control system
- 8. High pressure faucet drainage water
- 9. Over-temperature protection
- 10. Over-pressure protection
- II. Water shortage protection
- 12. Easy operate, safe and reliable



Specification

Model	DW-280D Digital	DW-280D Digital	DW-280D Digital
Volume	18L	24L	30L
Chamber size	Φ300*330mm	Φ300*460mm	Φ 300*520mm
Working temperature	115-129℃	115-129℃	115-129℃
Working pressure	0.145-0.165Mpa	0.145-0.165Mpa	0.145-0.165Mpa
Time scope	0-999min	0-999min	0-999min
Packing size	40*40*53cm	40*40*64cm	40*40*70cm
G.W/N.W	18kg/16kg	19kg/17kg	20kg/18kg
Power	AC220V/50Hz/2.0KW	AC220V/50Hz/2.0KW	AC220V/50Hz/2.0KW