

# DW-280A High Pressure Portable Autoclave

## Features

- 1. High Quality SUS304 stainless steel material
- 2. Electric heating system
- 3. Silicone seal
- 4. Dual scale indicate pressure gauge
- 5. Double safety valve to ensure safe
- 6. High pressure faucet
- 7. Easy operate, safe and reliable



## Specification

| Model               | DW-280A           | DW-280A           | DW-280A           |
|---------------------|-------------------|-------------------|-------------------|
| Volume              | 18L               | 24L               | 30L               |
| Chamber size        | Φ 300*330mm       | Φ300*460mm        | Φ 300*520mm       |
| Working temperature | 115-129℃          | 115-129℃          | 115-129℃          |
| Working pressure    | 0.145-0.165Mpa    | 0.145-0.165Mpa    | 0.145-0.165Mpa    |
| Packing size        | 40*40*53cm        | 40*40*64cm        | 40*40*70cm        |
| G.W/N.W             | 18kg/16kg         | 19kg/17kg         | 20kg/18kg         |
| Power               | AC220V/50Hz/2.0KW | AC220V/50Hz/2.0KW | AC220V/50Hz/2.0KW |