

## D Series Ultrasonic Cleaner

### Product Introduction

Ultrasonic cleaning machine is to use the cavitation effect of ultrasonic waves in the cleaning liquid to accelerate the dissolution and stripping of dirt, grease and other dirt in the inner hole and blind hole of the vessel being cleaned. DRAWELL ultrasonic cleaning machine adopts imported key equipment to ensure that the instrument has the advantages of fast cleaning speed, good effect, no damage to the surface and reduced labor intensity, etc.

### FEATURES

- 1.Digital display and setting ultrasonic cleaning time
- 2.Display working countdown time
- 3.Timing range: 1-600 minutes adjustable
- 4.Equipped with a special stainless steel basket, reduce sound cover
- 5.The shell and cover are made by stainless steel
- 6.The equipment controlled by Single-chip controller
- 7.The latest working parameters can be stored in memory after power-off



### Parameters

Model	Inner tank (mm)L/W/H	Volume (L)	Frequency (KHz)	Power (W)	Time (min)	Drainage
DW-50	145*80*65	0.8	40	50	1-120	/
DW-80	155*140*100	2	40	80	1-120	/
DW-100D	155*140*100	4.5	40	100	1-600	/
DW-120D	240*140*100	3	40	120	1-600	/
DW-120D	240*140*150	5	40	120	1-600	/
DW-3200D	300*155*150	6	40	180	1-999	/
DW-5200D	300*240*150	10	40	240	1-999	√
DW-5200D	300*240*150	10	40	360	1-999	√
DW-4200D	300*240*200	14.4	40	480	1-999	√
DW25-12D	500*300*150	22.5	40	600	1-999	√
DW25-12D	500*300*150	22.5	40	720	1-999	√
DW-800D	500*300*200	30	40	840	1-999	√