


MDJ Metallurgical Microscope



MDJ Metallurgical Microscope Size (mm)

Items	Specification	MDJ100	MDJ200
Eyepieces	WF 10×-18 mm	●	●●
Metallurgical LWD Infinity Plan Objectives	WF 10×-18 mm(adjustable,Reticule 0.1 mm)	○	○
	LPL 5×/0.15 W.D=15.00mm	●	●
	LPL 10×/0.25 W.D=18.48mm	●	●
	LPL 20×/0.40 W.D=8.35mm	○	○
	LPL 50×/0.75 W.D=2.55mm	●	●
	LPL 80×/0.80 W.D.= 0.85mm	○	○
	LWPlanFL 100×/0.85(Dry) W.D=3.0mm	○	○
Butterfly Binocular Head	Inclined 45°		●
Monocular Head	Inclined 45°	●	
Butterfly Trinocular Head	Inclined 45°, Built in 0.5X C-mount		○
Nosepiece	Quadplex	●	●
Mechanical Stage	Size: 180×180mm, Travel Range: 30×30mm, Round slide size: φ110mm	●	●
3W-LED Illumination	Input voltage : 100V~240V	●	●
12V/20W Illumination		○	○
Filters	Blue, Green, Grey	●	●
Micrometer	Division 0.01mm	○	○

Note: "●" In Table Is Standard outfits, "○" Is Optional Accessories.

Packing Size: 437mm×268mm×522mm Gross Weight: 8 kgs Net Weight: 7 kgs